



TESTIMONIAL

"Rodents will attack the bait as efficiently with a plastic wrap and the bait LASTS A LOT LONGER in the wrapper. From a cost analysis perspective for the pest controller, replacing bait regularly due to other pests cost money, but in this regard we find other pests go away and we only get hits from the intended rodent target only."

Charles Said-Vassallo

Commercial Pest Management



PRACTICALITY FEATURES

- MAKI® Block has a central hole so it can be easily secured inside a bait station
- High number of baits per bucket (375 for MAKI® Original, 275 for MAKI® Wrapped) allows more baiting points at times of low activity
- Convenient rectangular bucket for optimum storage

SAFETY FEATURES

- Bromadiolone is an anti-coagulant. Antidote Vitamin K1 is effective in case of accidental poisoning of children or pets
- MAKI Block contains Bifex, a bittering agent that reduces the risk of consumption by children or pets



MAKI® BLOCK

Available as



Liphatech

Product information contact:

Jerome Belime +61 421 946 857

email: info@liphatech.com.au

Disclaimer: Liphatech S.A.S. shall not be liable for any consequential or other loss or damage relating to the supply of or subsequent handling or use of this product, unless such liability is implied by statute or other law limited. All warranties, conditions or rights implied by statute or other law which may be lawfully excluded are so excluded. Where the liability of Liphatech S.A.S. for breach of any such statutory warranties and conditions cannot be lawfully excluded that may be limited to it re-supplying the product or an equivalent product or the cost of a product or an equivalent product, then the liability of Liphatech S.A.S. for any breach of such statutory warranty or condition is so limited.

® Registered trademark of Liphatech S.A.S

ALWAYS READ THE PRODUCT LABEL BEFORE USE

Available from:




MAKI® WRAPPED

WE'VE GOT YOU COVERED



OR



MAKI® ORIGINAL

IF YOU PREFER THEM UNDRRESSED (NAKED)

LIPHATECH

www.liphatech.com.au

LIPHATECH

www.liphatech.com.au



MAKI ORIGINAL

IF YOU PREFER THEM UNDRRESSED (NAKED)



YOUR EVERYDAY RODENT BAIT


- Contains the single-feed bromadiolone, originally developed by Liphatech.
- Highly mould-resistant, yet extremely effective.
- 20g extruded block. Use less than other bromadiolone baits.
- Ideal all-year round maintenance bait.



BEST USED:

- In commercial premises with ongoing rodent baiting programs to keep populations under control all year round
- Ideal for locations with high ambient humidity
- Can be used in a program to avoid resistance build-up in rodent populations.

- Registered on Black rats, Norway rats and House mice, for use in and around buildings
- Eg, alternate with a different bait containing a different active ingredient, such as Generation® Block.



HIGH MOULD-RESISTANCE

- Prolonged mould resistance for long-term efficacy

FIELD TRIAL



11 weeks, 15°C and 87% humidity

Tests run in a dairy factory with: 15 °C and 87% relative humidity (ie: baiting 12 Jan 2007 and assessment 03 April 2007)

MAKI WRAPPED

WE'VE GOT YOU COVERED



New concept!

BENEFITS:

- Plastic sachet increases longevity of baits protecting them from humidity and melting ==> less bait to replace
- Station stays cleaner, no mess at the bottom, other pests tend to not hibernate in the station

HOW TO USE:

- Mostly for commercial jobs, and for some residential situations (eg, subfloors, garden beds)
- Ideal when strong pressure from slugs/snails and/or insects (eg, cockroaches, etc).
- Very appropriate for audited sites; each sachet displays some regulatory information making Maki Block Wrapped the product of choice for professional pest managers
- During summer months, when internal bait station temperatures get really high

BEST USED:

- Protect your bait from other pests (eg slugs, snails, insects, etc).
- Highly weather-resistant. Less melt, less mould.
- No direct contact with bait for end-user.
- Information printed on sachet makes on-site data recording easier.

WOULD YOU LIKE IT WRAPPED?



Photo 1: bait preparation - Do not unwrap Maki Block.



Photo 2: bait placement inside station - station stays cleaner & consumption.



Photo 3: rodent gnawing on bait.



Photo 4: to secure bait, fold label edge on the block and punch it neatly on the top.



Photo 5: dead rat ==> bait consumed.



Photo 6: counter visit.

NOTE: Maki Block Wrapped will be performing best at times of low to medium rodent pressure. Replaces with an unwrapped block when rodent pressure increases.